

Quality Assurance in Radiodiagnostics Verification of Film Screen Contact |

according to EN ISO 4090



TEST DEVICE FFA 4090 R FOR RADIOGRAPHIC CASSETTES

Quality Assurance in Radiodiagnostics |

Verification of film screen contact and of the screen quality of radiographic cassettes/mammographic cassettes according to EN ISO 4090

A bad film-screen contact is causing a significant loss of sharpness of the x-ray. An important aspect of quality control – though often underestimated – therefore is the regular inspection of all cassettes in view of an even pressing of x-ray film and screen and also the screen integrity.

■ **This verification of the film-screen contact can be easily effected with:**

- the Test Device **FFA 4090 R** for radiographic cassettes
- the Test Device **FFA 4090 M** for mammographic cassettes

■ **Test procedure – short description – FFA 4090 R**

- Load the cassette equipped with intensifying screens with a new film (at least two minutes should elapse between cassette loading and exposure).
- Place the cassette on the table and adjust the x-ray field with the light-beam collimator, so that the cassette is fully illuminated.
- Place the test plate parallel to the film onto the cassette side which faces the X-ray tube.
- Set the following exposure parameters
Focus nominal value: max. 2 mm, Focus-to-screen distance: 150 cm,
X-ray tube voltage: max. 60 kV, Total filtration: 2 mm Al;
The exposure should be set so that within the open square measuring area, the optical density is 2.4 ± 0.4 .
- Release an exposure and develop the film.

■ **Test frequency**

By suspected changes in results or at least once annually

■ **Evaluation**

View the film on a film viewer from a distance of at least 2 m.
Collimate the film viewer to the film size.

■ **Result**

Good pressure is seen on the film by an even distribution of density over the whole film area. Spots with bad contact show a higher density and are seen as dark areas on the film. Spots with less density can be caused by damaged intensifying screens (e.g. due to chemical influences).