


Radcal

DAPcheck Plus

Quick and Easy Calibration Check - Of Installed DAP Meters

Verify X-ray DAP & Light Field Congruence



Traceable Measurements - The DAPcheck Plus is a reference class instrument for "field calibration" of patient dose measurement and control systems thus ensuring the validity of inter-institution patient dose comparisons.

FAST AND ACCURATE - Utilizing the high speed digitization of the Accu-Gold digitizer, the DAPcheck Plus displays DAP and DAP Rate of accumulated DAP on exposure completion. Designed as an accessory to the Accu-Gold.

DEPENDABLE - A tough ABS plastic housing protects the ion chambers and electronics that incorporate several patented features to ensure long term stability.

PLUS - Also provides X-ray to light field congruence



DAP VARIATION WITH KV			
Area (cm ²)	15 x 15	cm ² Dose rate	100
Time (ms)	200	DAP Rate (mSv/h)	200
Focus	100		
Field Size	100		

FIELD SIZE OF D20 FOR POC TESTING			
KV	Rad-1 DAP (µSv/cm ²)	Sample Size (cm ²)	% DAP VARIATION
10	1.2	1.0	±2%
15	1.8	1.0	±2%
20	2.4	1.0	±2%
25	3.0	1.0	±2%
30	3.6	1.0	±2%
35	4.2	1.0	±2%
40	4.8	1.0	±2%
45	5.4	1.0	±2%
50	6.0	1.0	±2%
55	6.6	1.0	±2%
60	7.2	1.0	±2%
65	7.8	1.0	±2%
70	8.4	1.0	±2%
75	9.0	1.0	±2%
80	9.6	1.0	±2%
85	10.2	1.0	±2%
90	10.8	1.0	±2%
95	11.4	1.0	±2%
100	12.0	1.0	±2%

THE GOLD STANDARD IN RADIATION MEASUREMENT

DAPcheck Plus KEY FEATURES AND BENEFITS:

KEY FEATURES

Complete DAP meter assessment:

Symmetrical Response:

Remote Control Software:

Optical and radiographic alignment markers:

X-ray to light field congruence:

BENEFITS

Measures DAP and DAP Rate over a full range of field sizes and beam qualities

Can be used with under couch tubes without the need for inversion

Automatic data capture with customizable templates

Setting reference field sizes made simple


Simultaneous measurements of light field congruence & DAP

XLGold (Excel Based) Customizable Template Example

OUTPUT CONSISTENCY			
kV	80	Baseline value Tolerance ± X%	
mA	200		
Time (ms)	100		
FCD (cm)	100		
Focus	Broad		
Sample	DAPCheckPlus Reading	Deviation from baseline	Press Ctrl-Shift-R to start Wait for cell to turn green before making exposure
1	1.163	0%	
2	1.164	0%	
3	1.164	0%	
4	1.162	0%	

DAP VARIATION WITH FIELD SIZE				
kV	80	X Company DAP units		µGy.m2
mA	200	DAP Tolerance ±X%		25%
Time (ms)	100			
Focus	Broad			
Press Cntl-Shift-R to start. Wait for cell to turn green before making exposure				
Area (cm)	X Company DAP (µGy.m2)	Sample	DAPCheckPlus (µGy.m ²)	% DAP VARIATION
10 x 10	10	1	10.56	-5%
15 x 15	25	2	25.21	-1%
20 x 20	46	3	46.40	-1%
30 x 30	105	4	105.30	0%


APPLICATIONS



DAP CALIBRATION

- DAP linearity with dose and field size.
- DAP calibration at the patient plane or at a reference distance.

X-RAY TO LIGHT FIELD CONGRUENCE SCALE



SPECIFICATIONS / TECHNICAL DATA:

Accuracy

DAP

Inclusive of all uncertainties (temperature, DAP pressure, rate, area and beam quality) ± 10%
Under reference conditions (10 mGy/min, 15 X 15 cm field, 80 kVp, 2.5 mm Al filtration) ± 7.5%

Digital Resolution

Dose area product rate

1µGy·m²/min to 0.91 Gy·m²/min

Rated range of use

Tube voltage

(40 – 150) kV

Automatic Temperature and Pressure Correction

Pressure

(80.0 – 106.0) kPa

Temperature

(+15 to +35) °C

Air humidity

(10 – 80) % rel. humidity (max. 20 g/m³)

Ionization chamber

Response versus radiation quality (50kV ... 150 kV, norm. to 100kV; acc. IEC 60580)

± 3 %

Quality equivalent filtration (70kV)

0.6 mm Al

Active area

max. (300 x 300) mm²

Dose area product

Connection to Accu-Gold Digitizer

Protection class (acc. IEC 60529)

One meter Disconnect Cable with 2 connectors
IP 41

Weight

2.32 kg; 5.11 lbs

Dimension

350 mm x 410 mm x 35 mm; 14 x 16 x 2 inches
(L x W x H)

All specifications subject to change.



PORTABLE CONVENIENT CARRY CASE

Foam Elevation Support Stand fits inside the interior of the Carry Case lid.

