

## Pro-Tec II<sup>®</sup> Syringe Shield

### Tungsten shielding

- Lead glass window, easy to replace
- Lightweight
- Fits most disposable syringes

The Pro-Tec<sup>®</sup> II Syringe Shield reduces hand exposure. It is designed for the clinician who has previously found it difficult to use conventional syringe shields. The Pro-Tec<sup>®</sup> II is thin, lightweight and easy to use. The barrel of the shield is constructed of .06" (1.5 mm) thick tungsten that will reduce radiation exposure from Tc-99m by 94%. A 5.05 density (1.5 mm lead equivalency) lead glass window provides protection and visibility. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content. A clip secures the lead glass to the tungsten barrel and makes it easy to replace scratched or broken glass. A stainless steel thumbscrew holds syringes firmly in place.

Pro-Tec<sup>®</sup> II Syringe Shields accommodate the standard sized 1 cc to 10 cc syringes.

Pro-Tec<sup>®</sup> II Syringe Shields with lead glass window:

**007-700** Syringe Shield, 1 cc

**007-800** Syringe Shield, 3 cc

**007-900** Syringe Shield, 5 cc

**007-909** Syringe Shield, 10 cc

**007-917** Syringe Shield, 20 cc



*Replacing scratched or broken glass is simple.*

MEDI-TEST SAS

Domaine Technologique - 4, rue René Razel - Immeuble Azur - F-91892 SACLAY CEDEX

Tel : +33 (0) 1 69 41 10 00 - Fax : +33 (0) 1 69 41 22 41 - [infos@meditest.fr](mailto:infos@meditest.fr)

Pro-Tec® II Syringe Shields without lead glass window:

**007-701** Syringe Shield, 1 cc

**007-801** Syringe Shield, 3 cc

**007-901** Syringe Shield, 5 cc

**007-911** Syringe Shield, 10 cc



*White coated interior makes it easy to read syringe contents.*

Replacement Glass for Pro-Tec® II Syringe Shields:

**127-700** Syringe Shield Replacement Glass, 1 cc  
For 007-700

**127-789** Syringe Shield Replacement Glass, 3 cc and  
5 cc  
For 007-800, -900, -835, -837, -903, -915

**127-999** Syringe Shield Replacement Glass, 10 cc  
For 007-909, -916

MEDI-TEST SAS

Domaine Technologique - 4, rue René Razel - Immeuble Azur - F-91892 SACLAY CEDEX

Tel : +33 (0) 1 69 41 10 00 - Fax : +33 (0) 1 69 41 22 41 - [infos@meditest.fr](mailto:infos@meditest.fr)