


GAFCHROMIC XR-CT dosimetry film 'Product Brief'



INTERNATIONAL SPECIALTY PRODUCTS



A BUSINESS UNIT OF ISP

Ideally suited to record density profiles associated with varying beam widths on CT scanners.


GAFCHROMIC XR-CT dosimetry film

Features



INTERNATIONAL SPECIALTY PRODUCTS

- Designed for the measurement of absorbed dose of low energy photons
- Substantial sensitivity enhancement with energy 20 kVp to 200 kVp.
- Dynamic dose range of 1 mGy to 500 mGy
- Energy response range from 80-140 kVp
- Self developing upon radiation
- Produces an instant density display
- Highly economical and reliable





A BUSINESS UNIT OF ISP

Structure of GAFCHROMIC® XR-CT dosimetry film



INTERNATIONAL SPECIALTY PRODUCTS





A BUSINESS UNIT OF ISP

GAFCHROMIC® XR-CT dosimetry film Sensitometric response



INTERNATIONAL SPECIALTY PRODUCTS



A BUSINESS UNIT OF ISP