

GAMMEX

MAMMOGRAPHY



Artifact Identification Phantom

GAMMEX 179

The Gammex 179 Artifact Identification Phantom is a full field device that allows for a grey film to be produced when the phantom is imaged. This will indicate any artifacts that may be caused by a mammographic system's grid or filters. It is also a great tool to verify that the film processor is working correctly. The Gammex 179 is a great QA check tool to verify that your mammographic system is working smoothly on all levels.

SPECIFICATIONS

Size. 3.8x30x24 cm (1.5x11.8x9.6 in)
Weight. 3.2 kg (7 lbs)



WWW.GAMMEX.COM

LASER ALIGNMENT

DIAGNOSTIC RADIOLOGY

ULTRASOUND

MAMMOGRAPHY