



Daily Laser and Light Field Plate

GAMMEX 443

The Gammex 443 Daily Laser and Light Field Plate is an effective tool for verifying the accuracy of the laser adjustments.

The Daily Laser and Light Field Plate determines that the lasers indicate the isocenter and the lines are straight, level, sharp, and in focus. The Gammex 443 provides a daily accuracy check of the following:

- Field size indicators
- Crosshair wander with collimator rotation
- Ceiling and wall laser alignment

SPECIFICATIONS

Size	20.3x20.3x3.6 cm (8x8x1.42 in)
Thickness	
Tool	3.18 cm (1.25 in)
Main Plate	0.64 cm (0.25 in) Side blocks
	extend 1.3 cm (0.5 in) above plate
Weight	0.31 kg (0.7 lbs)

