


NEW!



Beekley SPOTS[®] floral nipple markers

Beekley's new low-density, non-metallic nipple markers are perfect for digital mammography. The distinct 2.3mm round shape precisely and unmistakably documents the nipple's location on the image without artifact. Available in three popular nipple marker styles for your preference.



- image consistently with no artifact
- 2.3mm low-density, non-metallic pellets
- compatible with digital imaging equipment
- same quantity, same SPOTS[®] pricing
- distinct round shape permanently documents nipple's location
- soft, feminine floral designs for the less clinical, more patient-friendly environment



V-SPOTS[®]
 order code: 703
 contents per box: 130

- easy pick-up tab



N-FREE SPOTS[®]
 order code: 732
 contents per box: 110

- adhesive-free center



N-SPOTS[®] Soft 'n Stretchy
 order code: 790
 contents per box: 94

- pinch-free material and adhesive-free center

Try the FREE sample box enclosed!

Please check with your local distributor for pricing and product availability.

To locate a distributor in your area, visit www.beekley.com

Ask for Floral SPOT Markers only from Beekley!

BEEKLEY
CORPORATION

150 Dolphin Road, Bristol, CT 06010 Mailing: P.O. Box 369, Bristol, CT 06011 Tel: 1-800-233-5539 & +1-860-583-4700 Fax: 1-800-735-1234

0107