


Marking Specific Points and Lines of Reference on Treatment Sites

Conventional Simulation



Beekley Medical SPOTS® Metallic Markers: designed for use in conventional simulation and on port films. Lead-free metallic radiopaque pellet and linear markers are ideal for a variety of treatment planning points of reference.




T-SPOT .8mm wire used to identify size of abdominal tumor

- Lead-free, metallic pellet and linear markers image brightly
- A variety of different sizes in both pellet and linear shapes
- Identify and differentiate points of reference on treatment sites
 - > X-SPOT 1.5mm, Y-SPOT 2.0mm, V-SPOT 2.5mm pellets
 - > Z-SPOT 4.3mm burn-out resistant pellet
 - > T-SPOT .8mm continuous line
 - > S-SPOT .33mm continuous line
 - > D-SPOT radiopaque arrow

Perfect for Marking:


- ✓ Isocenters
- ✓ 3 point set-ups
- ✓ Outer canthus, nodes
- ✓ Scars
- ✓ Field borders, breast tangents
- ✓ Match lines

X-SPOT® Metallic Pellet




1.5mm pellets used to mark 3 point set-up on breast patient

Y-SPOT® Metallic Pellet



2.0mm pellets used to mark pelvic anterior and laterals for 3 point set-up on prostate patient

T-SPOT® and S-SPOT® Metallic Lines



.8mm and .33mm radiopaque continuous lines used to differentiate medial and lateral breast tangents

Beekley Medical SPOTS® Metallic Skin Markers

Designed for use in conventional simulation and on port films. Lead-free metallic radiopaque pellet and linear markers are ideal for a variety of treatment planning points of reference.



X-SPOT®

REF 101 QTY 150/box

- 1.5mm diameter pellet



Y-SPOT®

REF 102 QTY 132/box

- 2.0mm diameter pellet



V-SPOT®

REF 103 QTY 130/box

- 2.5mm diameter pellet



D-SPOT®

REF 104 QTY 99/box



Z-SPOT®

REF 114 QTY 40/box

- 4.3mm diameter pellet



T-SPOT®

REF 111 QTY 275cm cut to measure roll/box

- .8mm diameter line



S-SPOT®

REF 107 QTY 363cm cut to measure roll/box

- .33mm diameter radiopaque line

Product shown actual size

Call 1-800-233-5539 • Fax 1-800-735-1234
Visit www.beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1-860-583-4700 or email international@beekley.com

BEEKLEY MEDICAL™
when your diagnosis must be right™

One Prestige Lane, Bristol, CT 06010 Tel: 1-800-233-5539 or +1-860-583-4700 Fax: 1-800-735-1234 www.beekley.com

BEEKLEY, SPOT, SPOTS and X-SPOT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks.
BEEKLEY is a registered trademark in Canada.

© 2014 Beekley Corporation. All rights reserved.



REV: MS_SS_RT_0114