

GAFCHROMIC RTQA2 User Guide

GAFCHROMIC[®] RTQA2 is a self-developing, radiochromic film designed for QA and commissioning of radiotherapy modalities.

It is easy to use:

1. Position the film (with or without a phantom).
It does not matter which side of the film faces the source.
2. Expose the film,
3. See the image in real-time!

Position...



Position it, shoot, and see the exposure in real-time!

No processing needed.

...Shoot...



APPLICATIONS

GAFCHROMIC RTQA2 is currently available in four formats for:

- Light-field and radiation field alignment
- Precision star shots
- Position verification for HDR
- Autoradiography of implantable seeds, plaques, and other sources

...and Look!

RTQA2-1010
Light Field
Alignment Test

RTQA2-1010P
Pre-Scaled

RTQA2-1417
Large Field
Alignment Test



RTQA2-111
Brachytherapy Source Position Check

