

# LT-THERMOPLASTIC™

www.cdrrsys.ca  
403.850.7035

info@cdrrsys.ca  
1.855.856.7035


## LT Thermoplastic Immobilization & Sheets




### PRODUCT OVERVIEW

Our LT Thermoplastics are manufactured using the most current and advance processes available, ensuring the highest consistency with respect to stretch, working time and ultimate patient comfort.


All LT Thermoplastics are dry heat compatible. Ask your representative about our Dry Heat Convection Oven options.


**SL24-S-HYB-O**


Hybrid advanced open face mask. Compatible with all 3D surface monitoring systems


**SL32-SL 3.2 mm**


Drape and mold long LT thermoplastic. Available thicknesses: 2.4mm and 3.2mm


**SL32-O 3.2 mm**

Standard open face LT thermoplastic. Compatible with all 3D surface monitoring systems


**SL32-SXL 3.2 mm**


Drape and mold extra long LT thermoplastic, Available thicknesses: 2.4 mm and 3.2 mm


**SL24 2.4 mm**  
**SL32 3.2 mm**  
**SL48 4.8 mm**


Standard LT thermoplastic immobilization Available thicknesses: 2.4mm, 3.2mm and 4.8 mm. Compatible with shoulder tech


**STA32-01 3.2 mm**

Shoulder tech drape and mold LT thermoplastic



**SL24-L 2.4 mm**  
**SL32-L 3.2 mm**

Standard LT long thermoplastic immobilization. Available thicknesses in 2.4 mm and 3.2mm



**2.4 mm 3.2 mm**  
**RS24 (18"x24") RS32 (18"x24")**  
**RS24 (24"x36") RS32 (24"x36")**
**4.8 mm**  
**RS48 (17"x17")**  
**RS48 (12"x 24")**  
**RS48 (20"x 20")**

LT thermoplastic sheet, perforated or solid various sheet size. Available thicknesses: 2.4mm 3.2mm and 4.8 mm


Add **-01** to the order code for Solid sheets  
 Other sizes available, contact us for information


**SL24-XL 2.4 mm**  
**SL32-XL 3.2 mm**

Standard LT extra long thermoplastic immobilization. Available in thicknesses: 2.4mm and 3.2 mm


**SL24-S 2.4 mm**  
**SL32-S 3.2 mm**


Drape and mold head only LT thermoplastic. Available thicknesses: 2.4mm and 3.2mm. Compatible with shoulder tech


**2.4 mm 3.2 mm**  
**RS24P (12"x18") RS32P (12"x18")**  
**RS24P (18"x24") RS32P (18"x24")**  
**RS24P (24"x36") RS32P (24"x36")**

Pre-punched for handles, for various applications. Available thicknesses: 2.4mm, 3.2mm

Add **-01** to the order code for Solid sheets Other sizes available, contact us for information





CDR Systems

Precision Patient Positioning Products

[www.cdrsyst.ca](http://www.cdrsyst.ca) [info@cdsyst.ca](mailto:info@cdsyst.ca)

403.850.7035 1.855.856.7035

# HYBRID ADVANCED OPEN FACE MASK


The CDR Hybrid Advance Open Immobilization Mask is proven for all 3D surface monitoring systems.

Engineered for ultimate patient comfort, this mask never covers the patients face, making it especially ideal for claustrophobic patients.

Hybrid reinforcement layer maintains mask rigidity, compared to a full face mask

*Compatible with Intuition Frame™  
and MAYOMOLD™ Precision  
Immobilization*


*Compatible with Shoulder Tech  
LT Thermoplastic*



**SL24-S-HYB-O**



Hybrid Mask (SL24-S-HYB-O) shown with Shoulder Tech (STA-12) on Freedom System (LB-SBRT)



## WATER BATH




- ▶ LED digital microprocessor controlled water bath with precision temperature control within  $\pm 0.5^{\circ}\text{C}$  ( $\pm 1^{\circ}\text{F}$ )
- ▶ 120V/60hz, 1440 Watt heating element provides quick heat-up and even heat (North American specs)
- ▶ Met Lab approved to UL standards 60601-1 and CSA #C22.2NO.601
- ▶ Exterior Dim: 84cm x 66cm x 20cm (33"x26"x8")  
Interior Dim: 66cm x 56cm x 13cm (26"x22"x5")
- ▶ Side mounted drain with lockout for safe easy water exchange



**WB-LRG-2**  
Digital Water bath

**WB-LRG-2-EU**  
Digital Water bath  
(220V/50hz)



Water Drain



Micro Processor Control

## DRY THERMOPLASTIC OVEN



### CO-SM-1

Small compact convection Thermoplastic oven fits most LT Thermoplastic masks with manual control knobs and timer

*Exterior Dim: 12.5" x 20.75" x 15.5"*

*120V/60hz, 1750Watt, CE marked and NSF approved*



### CO-LRG-2

Larger convection Thermoplastic oven fits all LT Thermoplastic masks and ShoulderTech™. Thermoplastics digitally controlled

*Exterior Dim: 24" x 28" x 17"*

*Interior Dim: 19.5" x 15.5" x 10"*

*120V/60hz, 1450Watt, UL, NSF and CSA approved*

*Ask about our other thermoplastic dry heat ovens to fit your specific requirements*