


SHOULDER-TECH™ DRAPES FOR SHOULDER IMMOBILIZATION

The Shoulder-Tech™ Drape is designed to be used in conjunction with the Freedom™ System SRS and Head & Neck Module. Working in collaboration with the Intuition™ System and an LT-Thermoplastic Split Frame Drape type mask, the Shoulder-Tech™ Drape is easy to apply and provides shoulder immobilization while remaining extremely comfortable.




Shoulder Tech

<i>ST24-01</i>	LT-Thermoplastic Shoulder Tech 2.4mm	Pack of 5
<i>ST32-01</i>	LT-Thermoplastic Shoulder Tech 3.2mm	Pack of 5

FREEDOM PRONE HEAD & NECK MASK SET

Prone Mask set ensures correct support for prone patients. See page 23 for additional information on the Prone Immobilization System.



Freedom™ Prone Head and Neck Mask Set

<i>PH08-MSET</i>	LT-Thermoplastic Prone Head & Neck Mask set	Each
------------------	---	------

BL TYPE IMMOBILIZATION

BL Type LT-Thermoplastic immobilization masks are compatible with BrainLab™ type immobilization head only and head neck and shoulder immobilization platforms. Similar techniques and water bath temperature are implemented when working with this immobilization mask system.



BL Type Immobilization Masks

<i>BL406</i>	LT-Thermoplastic BL-Type A-Type	Pack of 5
<i>BL408</i>	LT-Thermoplastic BL-Type B-Type	Pack of 5
<i>BL410</i>	LT-Thermoplastic BL-Type Head, Neck & Shoulder Type	Pack of 5

MASK SHIMS

0.5mm Mask Shim for perfecting the mask fit. Supplied in packs of three and six CDR Mask Shims can be stacked to achieve customized correction.

Mask Shims

<i>MSHIM</i>	Face Mask Shims Thickness: 0.5mm	Pack of 6
<i>MSHIM-I</i>	Intuition Click™ Shims Thickness: 0.5mm	Pack of 3
<i>HSHIM</i>	Silverman Head Support Shims Thickness: 0.5mm	Pack of 3

