

**Multi-Modality Fiducial Markers:**

PinPoint® for Image Registration is a breakthrough fiducial marker with a unique conical design to assist in accurate image registration under both CT and MRI.



**PinPoint® for Image Registration**

order code: 128

contents per box: 30

Perfect for image registration performed under CT and MRI

Conical design with a 1.27mm diameter center hole

**MEDITEST SAS**

184 Rue Tabuteau - BP 80345 - 78533 BUC CEDEX  
Tel : +33 (0) 1 69 41 10 00 - Fax : +33 (0) 1 69 41 22 41  
infos@meditest.fr - www.meditest.fr